

nanoScan[®] MRI 3T and 7T

High-end MRI with the most robust cryogen-free superconducting magnet on the market





Founded 1990	Offices 7	Employees 300+	Publications 3200+	Countries 100+
------------------------	---------------------	--------------------------	------------------------------	--------------------------



Preclinical systems 400+	Clinical systems 1500+
------------------------------------	----------------------------------

About us


Mediso works in the field of medical imaging for 35 years with a profile of development, manufacturing, selling and servicing stand-alone and multi-modality imaging devices.

The company offers complete solutions from hardware design to evaluation and quantification software for clinical patient care and preclinical research.

Mediso has a leader position in the preclinical imaging market with over 400 commissioned systems around the world. Beyond the market leading nanoScan® PET/CT and SPECT/CT, Mediso also offers standalone MRI and integrated PET/MRI systems based on a cryogen-free magnet with 3T or 7T field strength and a PET insert for simultaneous PET/MRI imaging. Products are sold directly or through a distribution network in 100+ countries worldwide.



1990	1994	1998	2006	2009	2014	2016	2018	2019	2022	2022	2023
Mediso was founded	Introduction of the first Mediso gamma camera	Introduction of the Nucline™ X-Ring/4R, 4-detector dedicated brain SPECT	Launch of the first Mediso preclinical system, NanoSPECT/CT	Launch of preclinical nanoScan® PET/CT, world's first ever sub-mm resolution preclinical PET/CT	Introduction of MultiScan™ LFER 150 PET/CT, world's first sub-mm resolution mobile PET/CT	Introducing the nanoScan® PET/MRI 3T, world's first superconducting preclinical PET/MRI	Installation of the 100 th nanoScan® PET system	Introduction of nanoScan® PET/MRI with 7T MRI	Installation of the 300 th preclinical imaging system	Introduction of PET insert for nanoScan® PET/MRI 7T	Introduction of nanoScan® P123S PET/CT



Key features

3T AND 7T FIELD STRENGTH

100% Cryogen-free magnet

- No liquid helium or nitrogen
- Closed loop – no need to top-up helium

Wide-range of

- RF Coils
- Sequences

Compact design:

- **Small** footprint
- **Marginal** fringe field
- 480 / 970 kg (3T / 7T)

Powerful gradient

for DWI applications (up to 1050 mT/m)

Low-vibration, rear mounted

PulseTube Cryocooler for artefact free DWI-EPI

SmartMagnet™

- **Eco-friendly Idle Mode**
- **Active Quench Protection**



Upgrade possibility with 2-types of completely integrated PET systems

Designed for **dynamic studies**

- Freely accessible animal during the scan
- Minimized dead space for dynamic imaging
- Start dynamic acquisitions from touch screen (i.e. DCE),

Animal monitoring up to 3 animals

Easily accessible RF connection plate for surface coils with industry standard 50Ω connectors

User friendly hardware design


- Touchscreen can be used for tuning and matching
- Effortless RF coil removal
- PET Insert integrated with the RF coil

Easy to use, high-performance MRI platform

100% Cryogen-free magnet

The core of the nanoScan® MRI systems is the most robust **100% cryogen-free** superconducting magnet ever built for preclinical applications. It utilizes **conduction cooling** and **does not contain liquid helium or any other liquid cryogens** in any amount.


- It's base is a NbTi solenoid with multiple corresponding coils to maximize homogeneity and shielding thus reaching state-of-the-art homogeneity of ± 0.1 ppm @ 50 mm DSV and negligible fringe field outside the cryostat.
- Uniquely it features a **back mounted cryocooler** to significantly reduce conducted vibrations and to make maintenance easier.
- All electrically conductive cylindrical parts of the magnet were designed to minimize the residual eddy current after strong gradient pulses, this way achieving high-quality DWI images.



Unique back mounted cryocooler significantly reducing vibrations

SmartMagnet™ – Self-monitoring and management system

The patented* SmartMagnet™ technology enables one-click selection between different magnet modes.




*Discharge Controlled Superconducting Magnet and Superconducting Magnet Operating In Occasional Idling Mode


Vibration-free imaging

The system uses a **high-end pulse tube type cryocooler** that has significantly reduced mechanical vibrations compared to other cryocoolers due to the lack of moving mechanical components. This also results in

extended lifetime. The cryocooler is mounted on the back of the magnet and mechanically decoupled to further minimize vibrations inside the cryostat and therefore deliver ghost free EPI images.



2h long intensity stability measurement with single shot GRE-EPI




spinScan® next generation MRI spectrometer

The new Mediso spectrometer was **optimized for MRI applications** delivering an **ultra-low-noise** expandable **RF front-end** and **real-time dynamic shimming**.

The spectrometer is interfaced with three Mediso developed software to deliver a complete MRI workflow: the **Sequence Development Platform** for absolute control over the pulse sequences and reconstruction algorithms, the **FDA approved** and clinically validated **InterView™ FUSION** and **Nucline™** for routine acquisition planning, reconstruction, image post-processing and evaluation.

- **Dynamic shimming** for better EPI quality
- Eddy current compensation with **digital cross-term pre-emphasis**

- Highly expandable Rx and Tx channels
- Full **phased array** and **parallel Tx** support
- Continuous gradient monitoring
- Microsecond gradient resolution, for precise synchronization




Scalable Rx and Tx modules of the spinScan® spectrometer

High-end MRI applications made easy


Comprehensive pulse sequence library organized in application packages

The nanoScan® MRI systems are equipped with an extensive range of **continuously expanding, readily optimized protocols** for mice and rats including the most common MRI techniques as well as state-of-the-art pulse sequences and methods.


- **Anatomy:** Accurate anatomical localization and precise structural identification.
- **Angiography:** Voxel-wise Flow and Velocity analysis and dynamic flow measurement.
- **Cardiology:** High-fidelity Cardiac structural and functional mapping using trigger-based sequences.
- **CEST:** Saturation Transfer Analysis for Tissue Characterization.
- **Dynamic Imaging (fMRI and DCE):** Voxel-wise BOLD signal measurement for functional connectivity analysis. Quantitative perfusion and permeability Assessment with DCE.
- **Quantitative Mapping and Fat Water Imaging:** Precision Quantitative mapping & Fat-Water characterization.
- **UTE/ZTE:** Direct visualization of Rapidly Decaying Signals.
- **ASL:** Revealing perfusion abnormalities without contrast agent.
- Supports a wide range of **X-nuclei sequences** both imaging and spectroscopy.
- Improved SNR with **AI-based reconstruction**.
- **Diffusion:** Visualizing Pathological Alterations in Neurodegenerative- and vascular disorders.
- **Spectroscopy:** Quantitative biochemical analysis of tissues with metabolite ratio quantification.



Mouse, CINE Cardiac bright blood imaging



7T Rat brain, Single-shot SE EPI, 3 diffusion directions, BO (left) and corresponding ADC map (right)




In-vivo high resolution mouse brain imaging without (left) and with AI denoising reconstruction (right)

Achieving perfect SNR in every region with wide range of RF coils



Mediso offers a wide range of highly shielded, low noise RF coils designed to deliver the **best possible SNR**. The coils are **fully integrated** with the MultiCell™ system enabling precise animal positioning in relation to the coils, ensuring reproducible, quantitative results.

- **Transmit/receive volume coils for total body** imaging of up to obese rats or marmosets
- **Dedicated mouse and rat brain volume coils** with special imaging chambers
- **Flexible surface coils** of various diameters delivering excellent image quality and enhanced SNR
- **Phased array coils with multiple receiver channels** enabling parallel imaging for brain, heart or abdomen



Field upgradeable with best in-class PET systems


- **Full-scale PET-ring** on the front guaranteeing:
 - » unlimited quantitative imaging regardless the amount of radioactivity, acquisition time or animal size
 - » high-throughput by multiple-animal imaging capabilities, exploiting the large field-of-view
- **SiPM-based PET insert** ensures simultaneous PET/MRI studies without any interference between the PET and MRI readout.



MRI Applications

CINE Cardiac imaging at 7T


Improved spatial- and temporal resolution



ANIMAL MODEL: BALB/c mice
MRI ACQUISITION: CINE Cardiac Bright Blood, TR: 10.4 ms, TE: 1.9 ms
COIL: 2ch Cardiac surface coil

Rat Whole Body T1- and T2-weighted images at 7T


Rat Whole Body Imaging with T2 TSE (left) and T1 GRE (right) sequences



ANIMAL MODEL : Wistar rat
MRI ACQUISITION: Matrix size: 240 x 192, Voxel size: 0.4 x 0.4 x 1mm, NEX: 3
RF COIL: Rat Whole Body

Multi-Shot Spin Echo EPI at 7T


Multi-Shot (segmented) EPI imaging can achieve better resolution images by combining multiple EPI segments.



ANIMAL MODEL: BALB/c mouse
SEQUENCE: Multi-Shot EPI, TR:4000ms, TE: 35ms, FOV: 16mmx16mm, Matrix: 80x80, TH: 0.8mm
COILS: 36 mm Tx/Rx volume coil


Diffusion Tensor Imaging at 7T High-resolution ToF MRA at 3T

Single Shot DTI EPI in rat brain at 7T.



ANIMAL MODEL : Wistar rat
SEQUENCE: DW EPI 30 dir, b-value: 700, In-plane resolution: 230 µm, Slice Thickness: 1 mm, FOV: 23 mm x 23 mm x 23 mm, TR: 3750 msec; TE: 35 msec NEX: 2, Acquisition time: 13 min

With the nanoScan® MRI 3T and 7T exceptional quality ToF MRA images can be acquired without the need of any contrast agents. Even at 3T the ACA (Anterior Cerebral Artery), MCA (Middle Cerebral Artery) and CoW (Circle of Willis) is clearly visible.



ANIMAL MODEL : Wistar rat
SEQUENCE: 2D ToF MRA, In-plane resolution: 107µm, Slice Thickness: 100µm, FOV: 32 mm x 30 mm x 28 mm, TR: 70 msec, TE: 5.7 msec
COILS: 42mm Quadrature Rat brain coil and dedicated brain imaging chamber

Complete MRI workflow

Perform routine scans with the clinical validated Nucline™ acquisition software

Nucline acquisition software has been developed for **multimodal** medical imaging devices and is used in **clinical and preclinical** systems as well. It provides the same easy-to-use, integrated framework and main features for all the different modalities (**PET, SPECT, CT and MRI**). It integrates wide range of functionalities of acquisition, calibration, data management, reconstruction and visualization. Nucline has been developed with focus on **clean and user-friendly interface**, while complying to **industry standards** (CFR11, DICOM) and high level **cybersecurity** expectations.



1 PERSONALIZED ACCESS LEVELS


- **Routine:** A couple of clicks and the system is ready to run a study-specific, **optimized protocol**.
- **Advanced Tab:** Several acquisition and reconstruction parameters are editable providing the possibility of **further optimizing the protocols** for the study.
- **Research:** Access to all **system parameters** for researchers with significant experience



2 FOCUS ON QUALITY




- **Automatic, quick daily QC protocols**
- **Real-time diagnostic feedback**
- **Logged diagnostic data**



3 INTUITIVE GEOMETRIC PLANNER

- **Designing scans graphically** based on scout image or even any previous scan
- Setting up **advanced features** like shim box, saturation bands etc.
- **Real-time MRI signal** (selectable Real, Imaginary, Magnitude) during scan
- **Multi-Sequence monitor** to on-the-fly track progress of dynamic/gated/multi-series sequences
- **Easy-to-use image viewer** to quickly check the result image before next step



4 PREDEFINED, CONFIGURABLE PROTOCOLS


- Multimodality multi-step pre-saved **factory protocols** optimized for various applications
- Factory protocols can be copied, **edited, fine-tuned** by the User
- Study-optimized **User protocols** can be saved and loaded easily assuring **quick, reliable scanning**
- Protocols include **automated calibrations** (e.g. shim, RF, frequency etc.)
- Protocol steps can **run automatically** one by one
- Parameters are **validated** automatically

Steps	
Step	Status
(MRI) Shim3D Custom	✓
(MRI) SE EPI (4)	✓
(MRI) SE EPI (5)	✓
(MRI) SE EPI (6)	✓
(MRI) SE EPI (7)	✓
(MRI) DW EPI	✓
(MRI) DW EPI 2shot	✓


Analyze your quantitative data with the FDA approved InterView™ FUSION visualization and evaluation software

The FDA approved and clinically validated InterView™ FUSION multi-modal post-processing software is an essential part of system. It provides a wide range of functionalities to evaluate PET/SPECT/CT/MRI preclinical data for example:


- **Automatic MRI parametric evaluation**, e.g. T1, T2 and ADC map creation
- **AI-based MRI Denoising** reconstruction
- **Automatic multiple animal image separator**
- **Brain atlas**
- Wide range of 2D and 3D image viewers and rendering for visualization including 3D MIP and 3D Volume Rendering
- 3D and 4D data fusion via all image viewers and visualization of them over time frames
- Large variety of ROI/VOI tools
- **Time activity Curves (TAC)** of multiple ROIs/VOIs over 4D dynamic data with multiple statistics (min, max, mean, stdev, sum, etc.)




- **Automatic co-registration** procedures (rigid, affine and non-linear)
- **Advanced segmentation** methods
- Wide range of data input/output/export capabilities including video formats



Brain atlas



Automatic MRI parametric evaluation



Automatic multiple animal image separator



Complete MRI workflow

Have absolute control over your experiments with the **Sequence Development Platform**

The Sequence Development Platform enables you to **create your own sequences** or **modify the existing ones**.

With the intuitive GUI, programming and testing a new sequence is very easy and straightforward. A built-in tool helps you to test out various, freely customizable reconstruction and post-processing algorithms and adding new ones. Sequences can be directly run from this SW or exported to Nucline™ for routine scanning.


- Fully **interactive** Pulse-Sequence-Diagram
- RF and gradient pulse visualization with **unlimited custom shapes**
- Access to **factory sequence source codes**, GIT versioning
- Full access to raw data, K-space visualization
- Customizable **JAVA based code** for calculations and timing
- User defined **Python and MATLAB scripts** are executable as a reconstruction steps
- **Spectroscopy** module with advanced corrections
- Peak-picking and automatic curve fitting
- A variety of spectra visualization modes



Minimal installation requirements

The **compact** and **light-weighted** nanoScan® MRI systems can be **installed** and run practically in **any laboratory** due to their low installation and maintenance requirements.

As the systems are **100% cryogen-free** and **self-shielded** there is no need of any quench pipe or Faraday-cage. The 5-gauss line for both 3T and 7T models are close to the system, so they can be placed in really **small laboratories**, or other imaging systems like PET/CT or SPECT/CT can be easily placed in the same room. Moreover, in case of using the **silent air-cooled cryo-compressor** there is **no need for separate technical room**.



- **Light-weighted systems with small footprint**
480 kg / 970 kg (3T / 7T)
1050 / 2140 lbs (3T / 7T)
250 cm x 80 / 100 cm
- **Optional technical room** – In case of air-cooled cryo-compressor no separate technical room is needed
- **Post-processing workstation** – can be next to the acquisition workstation or at the researcher's room.
- **SiPM PET insert** – Optional upgrade for both 3T and 7T MRI systems.




nanoScan® MRI 7T reference installation with every system component (magnet, workstation, electronic cabinet, He-compressor, chiller) located in the same imaging room.

Animal handling

MultiCell™ imaging chambers


Mouse M

Inner space: 134×26 mm
Outer dimension: 463×32 mm
Up to 40 g



Mouse L (Standard)

Inner space: 141×31 mm
Outer dimension: 466×40 mm
Up to 80 g
Also available in BSL3 version



Rat L (Standard)

Inner space: 249×60 mm
Outer dimension: 580×70 mm
Up to 600 g



Mouse Triple


Inner space: 144×26 mm
Outer dimension: 488×70 mm
Up to 3×35-40 g



Monitoring and gating

- » ECG monitoring and triggering
- » Respiration monitoring and triggering
- » Temperature monitoring and control module
- » Accesible from touchscreen and workstation

Respiration and body temperature monitoring even up to four animals




Rat Dual

Inner space: 240×60 mm
Outer dimension: 590×70 mm
Up to 2×200 g



Mouse BSL-3

Inner space: 141×31 mm
Outer dimension: 578×60 mm
Up to 80 g



PrepaCell™

Supporting complete animal preparation before the scan, setting of:

- » Anaesthesia
- » Heating
- » Vital function monitoring

Eases workflow and increases throughput



nanoScan® MRI RF coils


Volume coil 34mm

Compatible with PET insert only
Available for 3T and 7T MRI systems
Inner diameter: 34mm





Volume coil 42mm

For mouse chambers
Available for 3T and 7T MRI systems
Inner diameter: 42 mm
X-nuclei versions also available



Volume coil 72 and 82 mm

For chambers available with MRI systems
Available for 3T and 7T MRI systems
Inner diameter: 42 mm
X-nuclei versions also available




Integrated brain and cardiac array coils

Compatible with 3T and 7T MRI systems
Compatible with mouse and rat MultiCell™ chambers
Available in 2 and 4 channel versions

Multifunctional flexible surface coils

Compatible with 3T and 7T MRI systems
Compatible with any MultiCell™ chambers
Available diameters: 10, 20, 30 mm



nanoScan® MRI 3T/7T
High-end MRI with the most robust cryogen-free magnet on the market



MRI

100% CRYOGEN-FREE • ROBUST MAGNET

- 3T and 7T field strength
- 100% Cryogen-free magnet
 - » No liquid helium or nitrogen
 - » Closed loop – no need to top-up helium
- Wide-range of RF Coils and Sequences
- Compact design:
 - » Small footprint
 - » Marginal fringe field
 - » 480 / 970 kg (3T / 7T)
 - » 1050 / 2140 lbs (3T / 7T)
- Powerful Gradient: (up to 1050 mT/m) for DWI application
- Low-vibration, rear mounted PulseTube cryocooler for artefact free DWI-EPI
- SmartMagnet™
 - » Eco-friendly idle mode
 - » Active quench protection
- Upgrade possibility with 2-types of completely integrated PET systems

nanoScan® PET/MRI 3T and 7T





Full-scale, quantitative PET combined with a robust, cryogen-free MRI



SPECT

HIGH SENSITIVITY • HIGH RESOLUTION • OUTSTANDING THROUGHPUT AT THE SAME TIME

- High resolution (0.3 mm *in vivo*) and high sensitivity 13 000 cps/MBq
- Largest field of view for large and multiple-animal imaging
- High-throughput
- Largest installation base 150+
- Highest flexibility:
 - » Wide isotope energy range, single or multiple: 20 keV – 1 MeV
 - » Various applications – optimized
- multi-pinhole collimators (e.g. MDP bonescan, dynamic, cardiac gated etc.)
 - » Animal models from tiny mouse up to large rabbit (6.5 kg)
 - » Different imaging schemes: helical, circular, full-stationary, 2D
- Parallel-hole collimators for imaging large animals
- List-mode acquisition





PET

BEST COUNT RATE PERFORMANCE • HIGHEST RESOLUTION WITH FREE ACCESS TO THE ANIMALS

- Highest resolution (< 0.7 mm)
- Largest transaxial field of view (12 cm)
- Longest axial field of view (up to 15 cm)
- Highest count rate performance (up to 1800 kcps @ 90 MBq / 2.43 mCi) supporting quantitative imaging in
 - » Radiotracer development
 - » Imaging of short half-life isotopes (e. g. ¹¹C, ¹³N, ¹⁵O)
 - » Multiple-animal imaging
- Free access to the animal supporting dynamic imaging
- Best NEMA sensitivity up to 10.5% (250-750 keV) and best Minimal Detectable Activity (MDA)
- Excellent quantification
- Largest installation base: 150+

nanoScan® SPECT/CT

Versatile SPECT/CT with absolute quantification and dedicated solutions for theranostic imaging



nanoScan® PET/CT

Real dynamic PET-system designed for quantitative studies



nanoScan® SPECT/CT/PET

Versatile SPECT with Real dynamic PET with absolute quantification



CT

HIGH POWER • HIGH RESOLUTION • LARGE FIELD OF VIEW

- High-resolution (30 µm)
 - Small voxel size (10 µm)
- Up to x7.6 zoom
- Variable transaxial field of view: 2–12 cm
- Highest power: 80 W X-ray tube for
 - » Large animals
 - » Better image quality
- Fast scanning
- Ex vivo samples
- Ultra-low dose protocol (<1 mGy for whole-body mouse)
- Real-time FBP and iterative reconstruction
- Respiration and cardiac gated reconstruction

MultiScan™ LFER 150 PET/CT

The ultimate tool for PET imaging in primates and medium sized animals



Specifications | nanoScan® MRI 3T and 7T

Magnet

Cryogen-free superconducting

Field strength

3T / 7T

Homogeneity

±0.1 ppm @ 50 mm DSV

Bore size

17 cm

Gradient coil inner diameter

101 mm

Gradient strength

Up to 1000 mT/m

Cryocooler

Back-mounted PulseTube

Quench protection

Yes, with SmartMagnet™

Faraday cage needed?

No, the system is self-shielded

Quench pipe needed?

No, the system is
100% cryogen-free

Rampable

Yes

400+ preclinical systems in
35 countries



nanoScan®
PET/CT

nanoScan®
SPECT/CT

nanoScan®
MRI 3T/7T

nanoScan®
PET/MRI 3T and 7T

nanoScan®
SPECT/CT/PET

MultiScan™
LFER150 PET/CT



Mediso Medical Imaging Systems
sales@mediso.com www.mediso.com

Headquarters
Budapest, Hungary

Global offices

USA and Canada

Arlington, VA
sales@medisousa.com

Belgium

Auderghem
info.belgium@mediso.com

United Kingdom and Ireland

Farnborough
info@mediso.uk

France

Strasbourg
info.fr@mediso.com

Germany and Austria

Münster
info@mediso.de

Poland

Łódź
biuro@mediso.pl



NS-M3T7T_0726_EN

Mediso reserves the right to change data without notice © Mediso 2026.

nanoScan® and spinScan® are registered trademarks of Mediso. MultiScan™, MultiCell™, Tera-Tomo™, PrepaCell™, Nucline™, InterView™, SmartMagnet™ are trademarks of Mediso.